C E 0086

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Toriimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

 EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
19 Loyang Way, #06-27 Singapore 508724
TEL: +65 6546 7588

PRINTED IN JAPAN

13007445
210231DK

GRADIA

- EN
- LIGHT-CURED COMPOSITE FOR CROWN & BRIDGE, INLAYS AND VENEERS
- DE
- LICHTHÄRTENDER KUNSTSTOFF ZUR HERSTELLUNG VON INLAYS, ONLAYS, KRONEN UND VENEERS, UND FÜR VERBLENDUNGEN VON KRONEN UND BRÜCKEN
- FR
- COMPOSITE DE LABORATOIRE PHOTOPOLYMERISABLE POUR COURONNES & BRIDGES, INLAYS ET FACETTES
- IT
- COMPOSITO FOTOPOLIMERIZZANTE PER CORONE, PONTI, INLAY E VENEER
- ES
- COMPOSITE FOTOPOLIMERIZABLE PARA CORONAS Y PUENTES, INCRUSTACIONES Y CARILLAS
- NL
- LICHTUITHARDEND COMPOSIT VOOR KRONEN, BRUGGEN, INLAYS EN VENEERS
- DA
- LYSPOLYMERISERENDE KOMPOSITMATERIALE TIL KRONER & BROER, INDLÆG OG FASADER
- SV
- LJUSPOLYMERISERANDE KOMPOSITMATERIAL FÖR KRON & BROTERAPI, INLÄGG OCH FASADER
- PT
- COMPÓSITO FOTOPOLIMERIZÁVEL PARA COROAS & PONTES, INLAYS e FACETAS
- EL
- ΦΩΤΟΠΟΛΥΜΕΡΙΖΟΜΕΝΗ ΣΥΝΘΕΤΗ ΡΗΤΙΝΗ ΓΙΑ ΣΤΕΦΑΝΕΣ & ΓΕΦΥΡΕΣ, ΕΝΘΕΤΑ και ΟΨΕΙΣ

For use only by a dental professional in the recommended indications.

INDICATIONS

- 1. Crown and Bridgework (with metal backing)
- 2. Metal free restorations: Jacket Crowns, Inlays, Onlays, Laminated Veneers
- 3. Implant superstructures

CONTRAINIDICATIONS

- 1. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- 2. In case of malocclusion, bruxism or clenching.

DIRECTIONS FOR USE

I. BUILDING UP PROCEDURE

- A. Anterior / Posterior Veneer Crown (with metal backing)
 - 1. Production of Metal Framework
 - a) Wax-up framework conventionally, occlusal surface should have a minimum of 1 mm thickness. Areas being covered by resin should have GC RETENTION BEADS II SSS added.
 - b) Sprue, invest, cast. Remove retainers in 0.5 mm width around the marginal veneer surface.
 - c) Clean and sandblast with 50 micron aluminum oxide.
 - d) Dry and clean with air gun. Do not clean with steam.
 - e) Immediately apply, with the brush, a coat of GC METALPRIMER II to retention areas. Allow to dry.
 - 2. Build up Veneer
 - a) Using the flat brush, immediately apply a layer of FOUNDATION OPAQUE (FO) over retention area (s). Light cure for 1 minute wit GC LABOLIGHT LV-III.
Note: For irradiation time by the GC visible light curing unit, refer to the chart given below.
 - b) Using the round brush, apply MARGIN OPAQUE (MO) in 1 mm width along the cervical line. Light cure for 1 minute.
 - c) Using the flat brush, apply two thin layers of OPAQUE (O) to the entire surface. Light cure each layer for 1 minute.
 - d) Progressively apply OPAQUS DENTIN (OD) as a cervical color in 2-3 mm width around the cervical area. Light cure for 3 seconds.
Note: In case of the crown with metal backing on the incisal edge, apply INTENSIVE COLOR (IC7, Lavender) along the incisal edge before applying OPAQUS DENTIN (OD). Light cure for 30 seconds.
 - e) Apply DENTIN (D), building up to the desired thickness. Light cure for 30 seconds. Use MAMELON STAIN (MS) or INTENSIVE COLOR (IC7) as required, and light cure for 30 seconds.
 - f) Apply ENAMEL (E) or FLOWABLE COMPOSITE (F-WT, F-WE) progressively from the incisal edge. Create a mamelon effect as required. Light cure for 30 seconds.
If desired, CERVICAL TRANSLUCENT (CT) can be applied to the cervical area, ENAMEL INTENSIVE (EI) to the incisal edge, HALO ENAMEL (HE) to the proximal surface. Light cure each layer for 30 seconds.
Note: FLOWABLE COMPOSITE should be covered by a layer of other pastes such as ENAMEL or ENAMEL INTENSIVE.
 - g) Coat the resin surface with GRADIA AIR BARRIER to ensure full polymerization of the surface and to avoid the oxygen inhibited layer.
 - h) Finally, light cure with the GC LABOLIGHT LV-III for 3 minutes. After light curing, remove the AIR BARRIER with water.
 - i) Grind & polish in usual manner. High-glaze finishing using GRADIA DIAPOLISHER on a felt wheel or leather buffing.
Note: For details of layering technique, refer to the GRADIA User's Manual.
- B. Anterior or Posterior Jacket Crown (metal free)
 - 1. Preparation of the working stone model
 - a) Prepare the stone die in usual manner.
 - b) Apply GRADIA DIE-HARDNER to the die stone surface.
 - c) Block out undercuts with wax.
 - d) Apply a layer of paraffin wax as die spacer as required.
 - e) Apply a thin layer of GRADIA SEPARATOR, a resin composite separator.
 - 2. Special pre-cure treatment
 - a) In case of a metal core, apply first a thin layer of MARGIN OPAQUE (MO) or OPAQUE (O). Light cure for 1 minute.
 - b) In case of a vital abutment tooth, apply first a thin layer of INTENSIVE COLOR (IC0). Light cure for 30 seconds.
 - 3. Build up of Crown
 - a) Apply a layer of OPAQUE (O) to the labial surface of the die and MARGIN OPAQUE (MO) to the lingual surface. Light cure each layer for 1 minute.
 - b) Apply SHOULDER DENTIN (SD) or DENTIN (one shade number just above the crown shade) in the cervical area. Light cure for 30 seconds.
 - c) For following build up on the labial surface, see instructions BUILDING UP PROCEDURE A. (Section 2. d-f) above.
 - d) For building up on the lingual surface, apply MAMELON STAIN (MS) to prevent the die from shining through the resin surface. Light cure for 30 seconds.
 - e) Apply OPAQUS DENTIN INTENSIVE (ODI) in the cervical area to cover 1/3 of the crown. Light cure for 30 seconds.
 - f) Progressively apply HALO ENAMEL (HE) from the incisal edge toward the cervical area. Also apply HALO ENAMEL (HE) to the proximal surface.
 - g) Continue BUILDING UP PROCEDURE A. (Section 2. g-i).
Note: For details of layering technique, refer to the GRADIA User's Manual.
- C. Posterior Inlay
 - 1. Preparation of the working stone model
Follow the instructions Preparation of the working stone model (B. Section 1. a-e) above.
 - 2. Special pre-cure treatment
Apply first a thin layer of INTENSIVE COLOR (IC0) to the cavity floor. Light cure for 30 seconds.
 - 3. Build up of Inlay
 - a) Build up CERVICAL TRANSLUCENT (CT) (highly translucent) as a dentin color in order to retain the natural tooth color. Light cure for 30 seconds.
Note: In the case of a stained tooth, it is recommended to first apply OPAQUE (O), MARGIN OPAQUE (MO) or OPAQUS DENTIN (OD) to the cavity floor. Light cure for O and MO for 1minute or OD for 30 seconds.
 - b) Build up ENAMEL INTENSIVE (EI) or ENAMEL (E), shape and contour. Light cure for 30 seconds.
 - c) Continue BUILDING UP PROCEDURE A. (Section 2. g-i) above.
Note: For details of layering technique, refer to the GRADIA User's Manual.

II. ADDITIONAL APPLICATION AND REPAIR

- 1. Roughen the resin surface using a carbide bur.
- 2. Apply a thin layer of GC COMPOSITE PRIMER to improve adhesion. Light cure for 1 minute with the GC LABOLIGHT LV-III.
- 3. Apply ENAMEL, TRANSLUCENT, etc. as required.

Prior to use, carefully read the instructions for use.

EN

- 4. Coat the added resin surface with GRADIA AIR BARRIER.
- 5. Light cure for 3 minutes.

SHADE COMBINATION CHART

Refer to the chart included in MASTER SET (in 10 shades) and STANDARD SET (in 6 shades).

IRRADIATION TIME

Light curing unit Paste	GC LABOLIGHT LV-III, II		GC STEPLIGHT SL-I	
	Pre-cure	Final cure	Pre-cure	Final cure
Foundation Opaque, Opaque, Margin Opaque	1 min.	—	—	—
Dentin, Opaqus Dentin, Opaqus Dentin Intensive, Shoulder Dentin, Enamel, Pearl Enamel, Halo Enamel, Enamel Intensive, Translucent, Cervical Translucent	30 sec.	3 min.	10 sec. *	—
Intensive Color, Mamelon Stain Flowable Composite	30 sec.	—	10 sec. *	—

* For one surface of a single crown

STORAGE

Store in a cool place away from direct sunlight. (Shelf life: 2 years from the date of manufacture)

PACKAGES

- 1. Set Packages
 - 1) MASTER SET (10 shades-A1, 2, 3, 3.5, B2, 3, 4, C2, 3, D3), 57 syringes with accessories
 - 2) STANDARD SET (6 shades-A2, 3, 3.5, B2, 3, C2), 26 syringes with accessories
 - 3) INTENSIVE COLOR SET (11 shades), 11 syringes
- 2. Refill Packages
 - 1) FOUNDATION OPAQUE, 2.4mL syringe (in 1 shade)
 - 2) MARGIN OPAQUE, 2.4mL syringe (in 1 shade)
 - 3) OPAQUE, 2.4mL syringe (in 22 shades)
 - 4) MAMELON STAIN, 2.4mL syringe (in 3 shades)
 - 5) INTENSIVE COLOR, 2.4mL syringe (in 15 shades)
 - 6) SHOULDER DENTIN, 2.9mL syringe (in 6 shades)
 - 7) OPAQUS DENTIN, 2.9mL syringe (in 23 shades)
 - 8) OPAQUS DENTIN INTENSIVE, 2.9mL syringe (in 6 shades)
 - 9) DENTIN, 2.9mL syringe (in 23 shades)
 - 10) ENAMEL, 2.9mL syringe (in 5 shades)
 - 11) HALO ENAMEL, 2.9mL syringe (in 1 shade)
 - 12) ENAMEL INTENSIVE, 2.9mL syringe (in 3 shades)
 - 13) PEARL ENAMEL, 2.9mL syringe (in 2 shades)
 - 14) TRANSLUCENT, 2.9mL syringe (in 5 shades)
 - 15) CERVICAL TRANSLUCENT, 2.9mL syringe (in 3 shades)
 - 16) FLOWABLE COMPOSITE, 2.4mL syringe (in 2 shades)
- 3. Accessories
 - 1) GRADIA SEPARATOR 5mL
 - 2) GRADIA DIE-HARDNER 5mL
 - 3) GRADIA DIAPOLISHER (solid) 8 g
 - 4) GRADIA AIR BARRIER 10mL
 - 5) GC COMPOSITE PRIMER 3mL
 - 6) GC METALPRIMER II 5mL
 - 7) SYRINGE PLUNGER
 - 8) SHADE GUIDE KIT (SG PLATE, SG CLIP FORMER)
 - 9) SG PLATE
 - 10) SG CLIP FORMER (12pieces)

CAUTION

- 1. In case of contact with oral tissues or skin, immediately remove with cotton. Flush with water.
- 2. In case of contact with eyes, flush immediately with water and seek medical attention.
- 3. GC METALPRIMER II, GC COMPOSITE PRIMER, GRADIA AIR BARRIER, GRADIA SEPARATOR and GRADIA DIE-HARDNER are flammable. Do not use near naked flame. Keep away from sources of ignition. Keep away from direct sunlight.
- 4. After dispensing, use immediately and replace syringe or bottle caps.
- 5. When adjusting or polishing the veneer surface, use a dust collector and wear a dust mask to avoid inhaling dust.

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This image shows a full page of white paper designed for handwriting practice. It features numerous horizontal dashed lines spaced evenly across the entire surface, providing a guide for letter height and placement. There are no margins, text, or other markings on the page.

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INDICATIONS

1. Crown and Bridgework (with metal backing)
2. Metal free restorations: Jacket Crowns, Inlays, Onlays, Laminated Veneers
3. Implant superstructures

CONTRAINIDICATIONS

1. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
2. In case of malocclusion, bruxism or clenching.

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1. Production of Metal Framework
 - a) Wax-up framework conventionally, occlusal surface should have a minimum of 1 mm thickness. Areas being covered by resin should have GC RETENTION BEADS II SSS added.
 - b) Sprue, invest, cast. Remove retainers in 0.5 mm width around the marginal veneer surface.
 - c) Clean and sandblast with 50 micron aluminum oxide.
 - d) Dry and clean with air gun. Do not clean with steam.
 - e) Immediately apply, with the brush, a coat of GC METALPRIMER II to retention areas. Allow to dry.
2. Build up Veneer
 - a) Using the flat brush, immediately apply a layer of FOUNDATION OPAQUE (FO) over retention area (s). Light cure for 1 minute wit GC LABOLIGHT LV-III.
Note: For irradiation time by the GC visible light curing unit, refer to the chart given below.
 - b) Using the round brush, apply MARGIN OPAQUE (MO) in 1 mm width along the cervical line. Light cure for 1 minute.
 - c) Using the flat brush, apply two thin layers of OPAQUE (O) to the entire surface. Light cure each layer for 1 minute.
 - d) Progressively apply OPAQUS DENTIN (OD) as a cervical color in 2-3 mm width around the cervical area. Light cure for 30 seconds.
Note: In case of the crown with metal backing on the incisal edge, apply INTENSIVE COLOR (IC7, Lavender) along the incisal edge before applying OPAQUS DENTIN (OD). Light cure for 30 seconds.
 - e) Apply DENTIN (D), building up to the desired thickness. Light cure for 30 seconds. Use MAMELON STAIN (MS) or INTENSIVE COLOR (IC7) as required, and light cure for 30 seconds.
 - f) Apply ENAMEL (E) or FLOWABLE COMPOSITE (F-WT, F-WE) progressively from the incisal edge. Create a mamelon effect as required.
Light cure for 30 seconds.
If desired, CERVICAL TRANSLUCENT (CT) can be applied to the cervical area, ENAMEL INTENSIVE (EI) to the incisal edge, HALO ENAMEL (HE) to the proximal surface. Light cure each layer for 30 seconds.
Note: FLOWABLE COMPOSITE should be covered by a layer of other pastes such as ENAMEL or ENAMEL INTENSIVE.
 - g) Coat the resin surface with GRADIA AIR BARRIER to ensure full polymerization of the surface and to avoid the oxygen inhibited layer.
 - h) Finally, light cure with the GC LABOLIGHT LV-III for 3 minutes. After light curing, remove the AIR BARRIER with water.
 - i) Grind & polish in usual manner. High-glaze finishing using GRADIA DIAPOLISHER on a felt wheel or leather buffing.
Note: For details of layering technique, refer to the GRADIA User's Manual.

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1. Preparation of the working stone model
 - a) Prepare the stone die in usual manner.
 - b) Apply GRADIA DIE-HARDNER to the die stone surface.
 - c) Block out undercuts with wax.
 - d) Apply a layer of paraffin wax as die spacer as required.
 - e) Apply a thin layer of GRADIA SEPARATOR, a resin composite separator.
2. Special pre-cure treatment
 - a) In case of a metal core, apply first a thin layer of MARGIN OPAQUE (MO) or OPAQUE (O). Light cure for 1 minute.
 - b) In case of a vital abutment tooth, apply first a thin layer of INTENSIVE COLOR (IC0). Light cure for 30 seconds.
3. Build up of Crown
 - a) Apply a layer of OPAQUE (O) to the labial surface of the die and MARGIN OPAQUE (MO) to the lingual surface. Light cure each layer for 1 minute.
 - b) Apply SHOULDER DENTIN (SD) or DENTIN (one shade number just above the crown shade) in the cervical area. Light cure for 30 seconds.
 - c) For following build up on the labial surface, see instructions BUILDING UP PROCEDURE A. (Section 2. d-f) above.
 - d) For building up on the lingual surface, apply MAMELON STAIN (MS) to prevent the die from shining through the resin surface.
Light cure for 30 seconds.
 - e) Apply OPAQUS DENTIN INTENSIVE (ODI) in the cervical area to cover 1/3 of the crown. Light cure for 30 seconds.
 - f) Progressively apply HALO ENAMEL (HE) from the incisal edge toward the cervical area. Also apply HALO ENAMEL (HE) to the proximal surface.
 - g) Continue BUILDING UP PROCEDURE A. (Section 2. g-i).
Note: For details of layering technique, refer to the GRADIA User's Manual.

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1. Preparation of the working stone model
Follow the instructions Preparation of the working stone model (B. Section 1. a-e) above.
2. Special pre-cure treatment
Apply first a thin layer of INTENSIVE COLOR (IC0) to the cavity floor. Light cure for 30 seconds.
3. Build up of Inlay
 - a) Build up CERVICAL TRANSLUCENT (CT) (highly translucent) as a dentin color in order to retain the natural tooth color. Light cure for 30 seconds.
Note: In the case of a stained tooth, it is recommended to first apply OPAQUE (O), MARGIN OPAQUE (MO) or OPAQUS DENTIN (OD) to the cavity floor. Light cure for O and MO for 1minute or OD for 30 seconds.
 - b) Build up ENAMEL INTENSIVE (EI) or ENAMEL (E), shape and contour. Light cure for 30 seconds.
 - c) Continue BUILDING UP PROCEDURE A. (Section 2. g-i) above.
Note: For details of layering technique, refer to the GRADIA User's Manual.

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1. Roughen the resin surface using a carbide bur.
2. Apply a thin layer of GC COMPOSITE PRIMER to improve adhesion. Light cure for 1 minute with the GC LABOLIGHT LV-III.
3. Apply ENAMEL, TRANSLUCENT, etc. as required.

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EN

4. Coat the added resin surface with GRADIA AIR BARRIER.
5. Light cure for 3 minutes.

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CAUTION

1. In case of contact with oral tissues or skin, immediately remove with cotton. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. GC METALPRIMER II, GC COMPOSITE PRIMER, GRADIA AIR BARRIER, GRADIA SEPARATOR and GRADIA DIE-HARDNER are flammable. Do not use near naked flame. Keep away from sources of ignition. Keep away from direct sunlight.
4. After dispensing, use immediately and replace syringe or bottle caps.
5. When adjusting or polishing the veneer surface, use a dust collector and wear a dust mask to avoid inhaling dust.

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