

Catégories	Normes
<p align="center"><u>Dents</u></p> SR Phonares II SR Vivodent S PE SR Orthotyp S PE SR Vivodent S DCL SR Orthotyp S DCL SR Vivodent SR Orthotyp Ivostar / Gnatostar	<p align="center">EN ISO 22112 Dentistry – Artificial teeth for dental prostheses</p>
<p align="center"><u>Résines</u></p> ProBase Cold / ProBase Hot SR Triplex Cold / Triplex Hot SR Ivocap SR Ivolen Ivobase High Impact / IvoBase Hybrid	<p align="center">EN ISO 20795-1 Dentistry – Base polymers – part 1: Denture base polymer</p>
<p align="center"><u>Résine provisoire</u></p> SR Ivocron	<p align="center">EN ISO 10477 Dentistry – Polymer-based crown and bridge materials</p>
<p align="center"><u>Composites</u></p> SR Adoro SR Nexco	<p align="center">EN ISO 10477 Dentistry – Polymer-based crown and bridge materials</p>
<p align="center"><u>Céramique sur alliage</u></p> IPS d.SIGN IPS InLine IPS Style Ceram	<p align="center">EN ISO 6872 Dentistry - Ceramic materials EN ISO 9693-1 Dentistry - Compatibility testing - Part 1: Metal-ceramic systems</p>
<p align="center"><u>Tout céramique</u></p> IPS e.max CAD IPS e.max Ceram IPS e.max Press IPS e.max ZirCAD IPS e.max ZirPress IPS Empress CAD IPS Empress Esthetic IPS Empress Universal	<p align="center">EN ISO 6872 Dentistry - Ceramic materials</p>
<p align="center"><u>Maquillants</u></p> IPS Ivocolor	<p align="center">EN ISO 6872 Dentistry - Ceramic materials</p>
<p align="center"><u>Alliages</u></p> 4all d.SIGN 30 Colado CC Colado NC	<p align="center">EN ISO 22674 Dentistry - Metallic materials for fixed and removable restorations and appliances</p>